

PURASORB® PDLG 7504

Product specification data sheet

Rev. No. 2 / June 2021

Description

PURASORB PDLG 7504 is a GMP grade copolymer of DL-lactide and Glycolide in a 75/25 molar ratio and with an inherent viscosity midpoint of 0.4 dl/g. It is supplied in the form of white to light tan granules. PURASORB PDLG 7504 is primarily used for drug delivery applications and is suitable for all commonly used formulation techniques.

Chemical composition75/25 DL-lactide/Glycolide copolymerMaterial numbers1001018327 / 1001018339 / (1840443)

 $\label{eq:molecular formula} \mbox{Molecular formula} \qquad \qquad [(C_6H_8O_4)_x(C_4H_4O_4)_y]_n$

Chemical name 3,6-dimethyl-1,4-dioxane-2,5-dione, polymer with 1,4-dioxane-2,5-dione

CAS Registry number 26780-50-7

Test	Method	Specification
Appearance	Visual test	White to light tan granules
Identity	FTIR spectroscopy	Conforms to reference
Monomer ratio, DL	NMR spectroscopy	72 - 78 mol %
Monomer ratio, G	NMR spectroscopy	22 - 28 mol %
Inherent viscosity	Viscometry Chloroform, 25 °C, c = 1 g/dl	0.35 - 0.45 dI/g
Specific rotation	Polarimetry Chloroform, 20 °C	(-2.0) - (+2.0) °
Water	Coulometric titration	max. 0.50 %
Mw	GPC, Chloroform	For information purpose only
Mn	GPC, Chlororform	For information purpose only
Mw/Mn	GPC, Chloroform	For information purpose only
Tin	ICP	max. 50 ppm
Residual solvent, total	GC, headspace	max. 0.10 %
Residual solvent, toluene	GC, headspace	max. 890 ppm
Residual monomer, total	GC	max. 2.0 %
Residual monomer, DL-lactide	GC	max. 2.0 %

corbion.com Page 1 of 2

Copyright © Corbion. All rights reserved. No part of this publication may be copied, downloaded, reproduced, stored in a retrieval system or transmitted in any form by any means, electronic, mechanical, photocopied recorded or otherwise, without permission of the publisher. No representation of warranty is made as to the truth of accuracy of any data, information or opinions contained herein or as to their suitability for any purpose, condition or application. None of the data, information or opinions contained herein may be relied upon for any purpose or reason. Corbion disclaims any liability, damages, losses or other consequences suffered or incurred in connection with the use of the data, information or opinions contained herein. In addition, nothing contained herein shall be construed as a recommendation to use any products in conflict with existing patents covering any material or its use.



Test		Method	Specification
Residual monomer, Glycolide		GC	max. 2.0 %
Elemental impurities		USP method 232	max. 10 ppm
Total microbial count		TAMC	max. 2 x 10 ³ CFU/g
Yeast and moulds		TYMC	max 2×10^2 CFU/g
Bacterial endotoxins		USP method 85	max. 6.00 EU/g
Acid number		Titration	max. 1.0 mg KOH/g
	out of four protective layers, an inner layer of clean room grade PE bag, an outer bag of aluminum coated polyester-PE laminate, inserted in the additional bag of PE for extra protection and shipped in sealed PE containers.		
Storage & Handling	When stored in the original packaging at low temperatures (-15°C), PURASORB PDLG 7504 keeps its initial properties for at least three years. (retest date) Stability studies indicate that, when stored in the original packaging at room temperature, PURASORB PDLG 7504 keeps its initial properties for at least one year.		
	Please refer to our current stability statement for most up to date information on storage stability.		
	Allow the material to reach room temperature before opening the packaging. After opening the original packaging PURASORB PDLG 7504 is best stored in an inert atmosphere and at low temperatures (-15°C).		